

## RIN3 Polyclonal Antibody

Catalog number: 12709-1-AP

Size: 20 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

RIN3; RIN3, Ras and Rab

interactor 3, FLJ22439,

FLJ11700, DKFZp762H1613



HepG2 cells were subjected to SDS PAGE followed by western blot with 12709-1-AP (RIN3 antibody) at dilution of 1:300

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Calculated RIN3 MW:	985aa, 108 kDa
Observed RIN3 MW:	108kd
Positive WB detected in:	HepG2 cells
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3390
GenBank accession number:	BC025248
Gene ID (NCBI):	79890
Full name:	Ras and Rab interactor 3

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.