

## DHH Polyclonal Antibody

Catalog number: 13889-1-AP


Size: 51 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

DHH; HHG-3; MGC35145



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 13889-1-AP(DHH Antibody) at dilution of 1:1000

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse skeletal muscle tissue
Calculated DHH MW:	396aa,44kd
Observed DHH MW:	42 kDa
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag4850
GenBank accession number:	BC033507
Gene ID (NCBI):	50846
Full name:	Desert hedgehog homolog (Drosophila)

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.