

APIP Polyclonal Antibody

Catalog number: 12089-2-AP


Size: 24 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

APIP; APAF1 interacting protein, APIP, APIP2, CGI 29, CGI29, dJ179L10.2, MMRP19, MTRu 1 P dehydratase



human kidney tissue were subjected to SDS PAGE followed by western blot with 12089-2-AP(APIP antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated APIP MW:	242aa, 27 kDa
Observed APIP MW:	27 kDa
Positive WB detected in:	Human kidney tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2729
GenBank accession number:	BC009077
Gene ID (NCBI):	51074
Full name:	APAF1 interacting protein

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.