

NGFRAP1 Polyclonal Antibody

Catalog number: 10446-1-AP

Size: 56 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

NGFRAP1; BEX3; Bex;

DXS6984E; HGR74; NADE



mouse liver tissue were subjected to SDS PAGE followed by western blot with 10446-1-AP(NGFRAP1 Antibody) at dilution of 1:400

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse liver tissue
Calculated NGFRAP1 MW:	13kd (22kd ,44kd)
Observed NGFRAP1 MW:	28kd
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP =

Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag0378
GenBank accession number:	BC003190
Gene ID (NCBI):	27018
Full name:	Nerve growth factor receptor (TNFRSF16) associated protein 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.