



GENERAL SPECIFICATIONS


| | |
|-----------------------|---------------------------|
| MATERIAL: | Polypropylene Homopolymer |
| COATING: | None |
| COLOR: | Natural |
| GRADUATION MARKS: | 250ul, 500ul, 750ul, 1mL |
| RACK CONFIGURATION: | Multichannel (8x12) |
| TEMPERATURE RANGE: | -80°C to 120°C |
| AUTOCLAVABLE: | 121°C at 15 psi, 15 min. |
| AUTOMATION COMPLIANT: | No |
| WORKING VOLUME: | 1.0mL |

CRITICAL DIMENSIONS

| | |
|-------------------|----------|
| LENGTH (OVERALL): | 100.00mm |
| COLLAR HEIGHT: | 28.50mm |
| BASE (OD): | 8.90mm |
| DISPENSING (ID): | 0.71mm |
| DISPENSING (OD): | 1.25mm |

PACKAGING CONFIGURATIONS

| | |
|--------------------------|--------------------------------------|
| PACKAGING QUANTITY (CS): | 96/Rack; 10 Racks/Unit; 5 Units/Case |
| CASE DIMENSIONS (IN.): | 19 x 10.75 x 24 |
| CASE WEIGHT (LBS.): | 30.5 |



FT1MXL-H96-NS

STORAGE

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone or ignition.

