

## GIPC2 Polyclonal Antibody

Catalog number: 12832-1-AP

Size: 92 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

GIPC2; FLJ20075, GIPC2,

SEMCAP 2, SEMCAP2



HEK-293 cells were subjected to SDS PAGE followed by western blot with 12832-1-AP(GIPC2 antibody) at dilution of 1:300

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated GIPC2 MW:	315aa,34 kDa
Observed GIPC2 MW:	34 kDa
Positive WB detected in	HEK-293 cells
Recommended dilution:	WB: 1:1000-1:10000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3566
GenBank accession number:	BC036075
Gene ID (NCBI):	54810
Full name:	GIPC PDZ domain containing family, member 2

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.