

AFF4 Polyclonal Antibody

Catalog number: 13197-1-AP


Size: 52 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

AFF4; AF4/FMR2 family member 4, AF4/FMR2 family, member 4, AF5Q31, AFF4, MCEF, Protein AF 5q31



Jurkat cells were subjected to SDS PAGE followed by western blot with 13197-1-AP(AFF4 Antibody) at dilution of 1:600

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated AFF4 MW:	365aa,40 kDa
Observed AFF4 MW:	150 kDa
Positive WB detected in	Jurkat cells, HL-60 cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3931
GenBank accession number:	BC025700
Gene ID (NCBI):	27125
Full name:	AF4/FMR2 family, member 4

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.