

High Quality, Clear Polypropylene Conical Bottom Centrifuge Tubes



- Leak Proof
- 95 kPa Compliant
- Temperature range from -80°C to 120°C
- Maximum RCF
 - 50mL tube = 15,500 x g
 - 15mL tube = 12,000 x g
- Polypropylene Tubes with HDPE Caps
- RNASE Free, DNASE Free and Pyrogen Free

Thomas Number	Product Number	Size (mL)	Quantity
21A00B685	BP1501	15mL Centrifuge Tube-Flat-Top-Rim Seal Cap-Sterile	20 racks of 25-500/case
21A00B686	BP1502	15mL Centrifuge Tube-Flat-Top-Rim Seal Cap-Sterile	20 bags of 25/pkg-500/case
21A00B687	BP5001	50mL Centrifuge Tube-Flat-Top-Rim Seal Cap-Sterile	20 racks of 25-500/case
21A00B688	BP5002	50mL Centrifuge Tube-Flat-Top-Rim Seal Cap-Sterile	20 bags of 25/pkg-500/case

