


Graphing



| Item No. | Qty | Description | Inventory Code | Material |
|----------|-----|--|----------------|----------|
| D1 | 1 | 2L SU Storage Bags | ACAX0012 | ULDPE |
| B1 | 2 | 1/4" 3/8" (6.4*9.6mm) L-Flex Tubing, Length 500mm | DBBX0063 | TPE |
| B2 | 1 | 1/4" 3/8" (6.4*9.6mm) Silicone Tubing (LePure), Length 200mm | DAAX0066 | Silicone |
| C1 | 1 | 1/4" (6.4mm) Hose Barb Non-Valved in-line coupling body | EBAX0005 | PC |
| C2 | 1 | MPC coupling plug | EBAX0028 | PC |
| C3 | 1 | 1/4" (6.4mm) Hose barb non-valved in-line coupling insert | EBAX0017 | PC |
| C4 | 1 | MPC coupling cap | EBAX0025 | PC |
| C5 | 1 | 1/4" (6.4mm) Female luer thread barb | EAXX0143 | PP |
| C6 | 1 | needle-less ampler | EABX0040 | PC |
| C7 | 1 | needle-less sampler cap | EABX0005 | PC |
| F1 | 6 | 2.5mm cable tie | ADAA0005 | Nylon |
| F2 | 3 | 3.6mm cable tie | ADAA0004 | Nylon |
| F3 | 3 | Robert clip (small) | EDBA0026 | PA |

| | | | | | | |
|--|--|---|--------------|--------------------|--------------|----------------|
| Customer Approval | | LePure Biotech Co., Ltd. | | | | |
| Company: | | Drawing No. HP002-4AK | | | | |
| Signature: | | Title: SU Storage Bags | | | | |
| Date: | | Designed By: | Date: | Checked By: | Date: | Edition |
| | | Jiajia Guo | 2023.03.12 | Alan | 2023.03.12 | Rev01 |
| <small>Design Approval the manufacturing of the system(s) based on this design. By signing, you also assume full</small> | | <small>We request that you do not disclose this information to any third party involved in this drawing, in order to protect our customer's confidentiality and intellectual property</small> | | | | |