

IMPROVACUTER® Gel and Clot Activator Tubes

IMPROVACUTER® Gel and Clot Activator Tubes are clinically proven with various clinical studies on a numerous analysis and ready to be used at the diagnostic labs. We value and focus on our patient's health and the accuracy of each test result.

Features:

1. Collection tube with spray-coated silica and polymer gel designed for effective serum separation.
2. Sterilization: Radiation, 10-6 SAL, ISO11137
3. Safety gold-color closure prevents contamination of serum from separated cellular components.
4. Blood clotting time: 30 minutes
5. Polyethylene terephthalate (PET) tube construction for reducing the risk of tube breakage during the shipping, handling, centrifuging at the high speeds between 1800-2200g.

Product code	Vol.(ml)	Tube Size(mm)	Additive	Color	Closure Type	Material	Pack/Carton Qty
623050202	5ml	13x100	Clot activator/Gel	Golden	Safety Cap	PET	100/1000

Application uses:

Venous Blood Collection tube.