

Hei-CHILL 1200 NR

Technical datasheet



Dimensions (W × H × D)	450 x 650 x 550 mm
Weight	N/A
Display	LED
Working temperature range	-10°C...40 °C
Display resolution	0.1°C
Interfaces	RS 232
Resolution	N/A
Pump flow max.	20 L/min
Noise Level	59 dB(A)
Pump pressure max.	1.3 bar
Power Supply	115 V; 60 Hz
Filling Volume (Min)	7 L
Filling Volume (Max)	14 L
Power Consumption (Max)	N/A
Power Plug	Cold Device Plug (IEC 60320, C14)
Protection Rating	IP 32
Temperature Stability	0.5 ± K
Refrigerant Stage 1	R-290 (GWP 3) A3; 0,085 kg; 0,0 t CO ₂ -eq
In/Outlet	G 3/4"
Safety Fittings	NFL

Resolution of Setting	0.1°C
Permissible Temperature Range During Storage	5°C...40°C
Permissible Temperature Range During Transport	-20°C...43°C
Protection Class	1