

ARPP-21 Polyclonal Antibody

Catalog number: 11829-1-AP

Size: 57 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

ARPP-21; ARPP 21, ARPP21,

ARPP-21, FLJ32997, RCS,

TARPP



mouse brain tissue were subjected to SDS PAGE followed by western blot with 11829-1-AP(ARPP-21 antibody) at dilution of 1:800

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated ARPP-21 MW:	89kd
Observed ARPP-21 MW:	86 kDa
Positive WB detected in	Mouse brain tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2399
GenBank accession number:	BC017805
Gene ID (NCBI):	10777
Full name:	Cyclic AMP-regulated phosphoprotein, 21 kD

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.