

TULP2 Polyclonal Antibody

Catalog number: 12668-1-AP

Size: 24 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

TULP2; Cancer/testis antigen
65, CT65, tubby like protein 2,

Tubby related protein 2,

TUBL2, TULP2



HeLa cells were subjected to SDS
PAGE followed by western blot with
12668-1-AP(TULP2 antibody) at
dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Calculated TULP2 MW:	520aa, 59 kDa
Observed TULP2 MW:	59kd
Positive WB detected in	HeLa cells, HEK-293 cells, human liver tissue, human lung tissue, Jurkat cells
Recommended dilution:	WB: N/A

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3350
GenBank accession number:	BC026070
Gene ID (NCBI):	7288
Full name:	Tubby like protein 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.