

Technical specification for blackKnights 1000 µl

article no.: 49003-0000 = 10 trays with 96 tips in cardboard packaging
49003-1002 = 10 trays with 96 tips in cardboard packaging, sterilized
49003-1102 = 10 trays with 96 filter tips, in cardboard packaging, sterilized
49003-1104 = 5 trays with 96 filter tips in blister packaging, sterilized

aliquot tips: 96 disposable tips in tray (SBS format)

level sensing: yes (resistance < 30 kΩ)

compatible: Tecan Genesis, Freedom Evo

colour: black (tray = red)

aliquot tip size: 1000 µl

aliquot volume: 50 – 1000 µl

length: 95,5 mm

coaxiality: < 2 mm

outside diameter: 7.8 mm

CE: IVD-compliant according 98/79EC, annex II

purity tested: Dnase /Rnase-free, endotoxin-free

Made in Germany under clean room conditions