

## Synaptotagmin-9 Polyclonal Antibody

Catalog number: 12793-1-AP

Size: 45 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

SYT9; FLJ45896, Synaptotagmin

9, synaptotagmin IX, SYT9,

SytIX



human brain tissue were subjected to SDS PAGE followed by western blot with 12793-1-AP(SYT9 antibody) at dilution of 1:500

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated Synaptotagmin-9 MW:	491aa,52kd
Observed Synaptotagmin-9 MW:	52kd
Positive WB detected in	Human brain tissue, human kidney tissue, mouse eye tissue
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3520
GenBank accession number:	BC029605
Gene ID (NCBI):	143425
Full name:	Synaptotagmin IX

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.