

FGFR2 Polyclonal Antibody

Catalog number: 13042-1-AP

Size: 24 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

FGFR2; BEK, BFR 1, CD332,


CEK3, CFD1, ECT1, FGFR 2,

FGFR2, FLJ98662, JWS, K SAM,

KGFR, KSAM, TK14, TK25



A549 cells were subjected to SDS PAGE followed by western blot with 13042-1-AP(FGFR2 antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using 13042-1-AP(FGFR2 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Applications

Tested applications:	ELISA, WB, IF
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated FGFR2 MW:	709aa,79 kDa
Observed FGFR2 MW:	120 kDa,85-92 kDa
Positive WB detected in	A549 cells, A431 cells, HeLa cells, HepG2 cells, SGC-7901 cells
Positive IF detected in	HepG2 cells
Recommended dilution:	WB: 1:200-1:2000 IF: 1:10-1:100

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3701
GenBank accession number:	BC039243
Gene ID (NCBI):	2263
Full name:	Fibroblast growth factor receptor 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.