

## NPY Polyclonal Antibody

Catalog number: 12833-1-AP


Size: 30 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

NPY; PYY4



Immunohistochemistry of paraffin-embedded human pancreas cancer using 12833-1-AP (NPY Antibody) at Dilution 1:50 (under 40x lens)

### Applications

Tested applications:	ELISA, IHC
Cited applications:	IHC
Species specificity:	Human, mouse, rat; other species not tested.
Cited species:	Human
Calculated NPY MW:	97aa,11kd
Observed NPY MW:	
Positive IHC detected in:	Human pancreas cancer
Recommended dilution:	IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag3568
GenBank accession number:	BC029497
Gene ID (NCBI):	4852
Full name:	Neuropeptide Y

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.