

ULBP2 Polyclonal Antibody

Catalog number: 13133-1-AP

Size: 22 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

ULBP2; ALCAN alpha, N2DL 2,

N2DL2, NKG2D ligand 2,

NKG2DL2, RAET1H, UL16

binding protein 2, ULBP2



Jurkat cells were subjected to SDS PAGE followed by western blot with 13133-1-AP(ULBP2 antibody) at dilution of 1:400

Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Calculated ULBP2 MW:	246aa,27 kDa
Observed ULBP2 MW:	30-35 kDa
Positive WB detected in	Jurkat cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3798
GenBank accession number:	BC034689
Gene ID (NCBI):	80328
Full name:	UL16 binding protein 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.