

Universal Pipet Tip Compatibility

Following are some of the pipet brands and series generally compatible with ILP Filtered Universal Pipet Tips.

BIOHIT	PROLINE Single channel PROLINE Multi channel E-LINE Single channel E-LINE Multi channel M-LINE Multi channel	BRANDTECH	TRANSFERPETTE Single channel TRANSFERPETTE Multi channel
CAPP AERO	Pipettor Single channel Pipettor Multi channel	EPPENDORF	ELECTRONIC REFERENCE RESEARCH Single channel RESEARCH Multi channel
THERMOLABSYSTEMS	FINNPIPETTE Single channel	GILSON	PIPETMAN Single channel
NICHIRYO/OXFORD	Benchmate	RAININ	EDP3 Single-channel EDP3 Multi-channel LATCH MODE
SOCOREX	CALIBRA/ACURA Single channel CALIBRA/ACURA Multi channel	Thomas Scientific	SIGNATURE Single channel