

PIGA Polyclonal Antibody

Catalog number: 13679-1-AP


Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

PIGA; GPI3; PIG-A



HeLa cells were subjected to SDS PAGE followed by western blot with 13679-1-AP(PIGA Antibody) at dilution of 1:1000

Applications

Tested applications: ELISA, WB
Species specificity: Human, mouse, rat; other species not tested.
Positive WB detected in: HeLa cells, HEK-293 cells, Jurkat cells
Calculated PIGA MW: 484aa,54kd
Observed PIGA MW: 38kd
Recommended dilution: WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen: Ag4530
GenBank accession number: BC038236
Gene ID (NCBI): 5277
Full name: Phosphatidylinositol glycan anchor biosynthesis, class A

Product information

Purification method: Antigen affinity purification
Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.