

DNAL4 Polyclonal Antibody

Catalog number: 10388-1-AP


Size: 37 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

DNAL4; PIG27



mouse testis tissue were subjected to SDS PAGE followed by western blot with 10388-1-AP(DNAL4 Antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse testis tissue, HEK-293 cells
Calculated DNAL4 MW:	12kd
Observed DNAL4 MW:	12kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP =

Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag0587
GenBank accession number:	BC002968
Gene ID (NCBI):	10126
Full name:	Dynein, axonemal, light chain 4

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.