

Product Name:	ABS LD 50
Description:	ABS LD Class 50 Liter Liquid Nitrogen Capacity 0.49 Liter/Day Static Evaporation Rate Neck Diameter (mm): 63.5 Exterior Height (mm): 774.4 Exterior Diameter: 431.8 Designed for efficient storage of LN2 Low LN2 consumption Durable composite lid High strength neck tube reduces LN2 loss Strong, lightweight aluminum construction Safer operation & maintenance One-year parts & labor warranty Three-year vacuum warranty

