



**NEW!**

# DispensMate Plus

Bottle-Top Dispensers



Quick and easy volume adjustment

Smooth piston design provides effortless dispensing

Closure cap protects against reagent contact and is easy to attach and remove, even while wearing gloves.



Adapters allow fitting to standard reagent bottles



Filling tube adjusts easily to different size bottles

## Ordering information

### DispensMate Plus

Thomas No.	Cat No.	Volume ml	Graduation ml	A ≤ ±		CV ≤	
				%	ul	%	ul
1226W01	DZ7440505	0.5-5.0	0.10	0.5	25	0.1	5
1226W02	DZ7440110	1.0-10.0	0.20	0.5	50	0.1	10
1226W03	DZ7442525	2.5-25.0	0.50	0.5	125	0.1	25
1226W04	DZ7440550	5.0-50.0	1.00	0.5	250	0.1	50

### Features:

- Excellent chemical resistance
- Entirely autoclavable at 121°C
- Four ranges of bottle top dispenser covers a volume range from 0.5ml to 50ml
- Easy to clean and maintain
- The optional flexible discharge tube with safety handle permits fast and precise dispensing
- Made of PTFE, ETFE, FEP, BSG, PP
- 2.2g/cm<sup>3</sup> vapor pressure Max. 500mbar, viscosity max. 500mm<sup>2</sup>/s, temperature max. 40°C, density max. 2.2g/cm<sup>3</sup>
- Dispensmate's are supplied with S40, GL32, GL38, GL25, GL28 sized adapters.

**DRAGONLAB**