




General Purpose Sealing Film



Sealing Film, General Purpose
For sealing both polystyrene and polypropylene plates

Application:
Ideal for ELISA, Incubation and General Use

Temperature Range:
-40°C to +120°C

Catalog #	Description	Quantity
GP-PP-100	General Purpose Film	100 Sheets
GPS-PP-100	General Purpose Film, Sterile	100 Sheets

Features and Benefits

- ▶ Certified DNase- & RNase- free
- ▶ Designed to be used with both polystyrene and polypropylene plates
- ▶ 2.0mil clear polypropylene film features good clarity and easy handling
- ▶ Specially formulated adhesive provides an excellent seal to minimize evaporation, reduce cross-contamination and prevent spillage
- ▶ Removes easily and cleanly from plate surfaces
- ▶ Minimizes aerosols upon film removal
- ▶ Features two end-tabs with perforations for ease of handling, application and removal