

Specification Sheet

PRODUCT: 100mM Tris Hydrochloride, pH 8.2

CATALOG NO.: 472910 and
492910

LOT NO.: TBD

MANUFACTURE DATE TBD

EXPIRATION DATE: 1 Year from the date
of manufacture

STORAGE CONDITIONS: 4 - 40°C

Test Results

Test	Specification	Test Method	Result
Appearance	Clear/Colorless Liquid	Visual	TBD
pH	8.0 - 8.4	Potentiometric	TBD
Protease	None Detected	Enzymatic	TBD
DNase	None Detected	Fluorometric	TBD
RNase	None Detected	Fluorometric	TBD
Peroxides Content	≤ 5ppb	Enzymatic	TBD
Sterility	No Growth	USP <71>	TBD
Endotoxins	≤ 1 EU/mL	USP <85>	TBD

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

Manufacturing Site

Lifecycle Biotechnologies, 5858 Wright Drive, Loveland, Colorado 80538

Simple Lifecycle