



Gelatin Phosphate Buffer

M1359

Gelatin Phosphate buffer is recommended for toxin detection in food products when *Clostridium botulinum* is suspected.

Composition**

Ingredients	Gms / Litre
Sodium dihydrogen phosphate	4.000
Gelatin	2.000
Final pH (at 25°C)	6.2±0.1

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 6.0 grams in 1000 ml distilled water. Boil to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and dispense as desired.

Principle And Interpretation

Botulinum toxin (botox) types A-G are produced by heterogeneous strains of *Clostridium botulinum*. Botox types A,B,E and F have caused serious and sometimes fatal, cases of food borne illness in humans. The vast majority of botulinum outbreaks in red meat and poultry products have involved either toxin A or B.

The current botulinum toxin test method is the mouse bioassay procedure (1). Gelatin Phosphate Buffer is one of the reagent used in this test method.

Quality Control

Appearance

Cream to yellow coloured homogeneous coarse powder

Gelling

Semisolid comparable with 0.2% gelatine

Colour and Clarity of prepared medium

Colourless clear solution forms in tubes.

Reaction

Reaction of 0.6% w/v aqueous solution at 25°C. pH : 6.2±0.1

pH

6.10-6.30

Cultural Response

M1359: Cultural characteristics observed after an incubation at 35-37°C for 24-48 hours.

Organism	Inoculum (CFU)	Growth	Gelatinase reaction
Cultural Response <i>Clostridium botulinum</i> ATCC 25763	50-100	good	Positive reaction

Storage and Shelf Life

Store below 30°C in tightly closed container and use freshly prepared medium. Use before expiry date on the label.

Reference

1.L. Victor Cook, Wei Hwa, Charles P. Lattuade and Gerri M. Ransom Methods for the detection of *Clostridium botulinum* toxins in meat and poultry products. USDA/FSIS Microbiology Laboratory Guidebook, 3rd edition, Chapter 14, 1998.

Revision : 1 / 2011



Technical Data



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.