

SLITRK4 Polyclonal Antibody

Catalog number: 13835-1-AP


Size: 29 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

SLITRK4; DKFZp547M2010



mouse thymus tissue were subjected to SDS PAGE followed by western blot with 13835-1-AP(SLITRK4 Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse thymus tissue
Calculated SLITRK4 MW:	837aa,94kd
Observed SLITRK4 MW:	96kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag4781
GenBank accession number:	BC040986
Gene ID (NCBI):	139065
Full name:	SLIT and NTRK-like family, member 4

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.