

## BMP5 Polyclonal Antibody

Catalog number: 13253-1-AP


Size: 34 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

BMP5; BMP 5, BMP5, bone morphogenetic protein 5



mouse lung tissue were subjected to SDS PAGE followed by western blot with 13253-1-AP( BMP5 Antibody) at dilution of 1:600

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse; other species not tested.
Calculated BMP5 MW:	454aa,52 kDa
Observed BMP5 MW:	52 kDa
Positive WB detected in	Mouse lung tissue, HUVEC cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3952
GenBank accession number:	BC027958
Gene ID (NCBI):	653
Full name:	Bone morphogenetic protein 5

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.