

SAPK4 Polyclonal Antibody

Catalog number: 10217-1-AP

Size: 22 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

MAPK13; MAP kinase 13, MAP

kinase p38 delta, MAPK 13,

MAPK13, p38delta, PRKM13,

SAPK4



HepG2 cells were subjected to SDS PAGE followed by western blot with 10217-1-AP (MAPK13 Antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated SAPK4 MW:	38 kDa
Observed SAPK4 MW:	42 kDa
Positive WB detected in:	HepG2 cells
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag0255
GenBank accession number:	BC000433
Gene ID (NCBI):	5603
Full name:	Mitogen-activated protein kinase 13

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.