



Flip-Up Disposable Face Shields



This disposable face shield provides excellent visibility while protecting the face and neck from hazardous splashes and spray.

- Optically clear for peripheral visibility. One size fits most.
- Limit the spread of viruses, bodily fluids, and health hazards
- 1" wide plastic headband strap
- Adjustable rear hook & loop fastener fits head circumference 17" - 26"
- 14"W x 8"H x .020" thick plastic front
- Plastic face shield is non-sterile
- Face shield features flip-up and down with side aluminum rivets

Thomas No.	Description	Qty
20A00H653	Face Shield, Flip-Up Disposable, Clear	50/Cs

For more information, visit [ThomasSci.com](https://www.ThomasSci.com) or contact your Thomas Scientific Sales Representative today!

Disclaimer: Pricing and offers subject to change without notice.