

## CDC20B Polyclonal Antibody

Catalog number: 13376-1-AP

Size: 28 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

CDC20B; CDC20B, FLJ37927,

G6VTS76519

### Applications

**Tested applications:** ELISA  
**Species specificity:** Human; other species not tested.  
**Calculated CDC20B MW:** 192aa,22 kDa  
**Observed CDC20B MW:**

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

**Immunogen:** Ag4198  
**GenBank accession number:** BC037547  
**Gene ID (NCBI):** 166979  
**Full name:** Cell division cycle 20 homolog B (*S. cerevisiae*)

### Product information

**Purification method:** Antigen affinity purification  
**Storage:** PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.