

MEOX2 Polyclonal Antibody

Catalog number: 12449-1-AP


Size: 69 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

MEOX2; GAX, Homeobox protein MOX 2, MEOX2, mesenchyme homeobox 2, MOX2



IP Result of anti-MEOX2 (IP:12449-1-AP, 4µg; Detection:12449-1-AP 1:300) with human placenta tissue lysate 1520µg.



human placenta tissue were subjected to SDS PAGE followed by western blot with 12449-1-AP (MEOX2 Antibody) at dilution of 1:600

Applications

Tested applications:	ELISA, IP, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated MEOX2 MW:	303aa, 33 kDa
Observed MEOX2 MW:	33 kDa
Positive WB detected in:	Human placenta tissue
Positive IP detected in:	Human placenta tissue
Recommended dilution:	WB: 1:200-1:2000 IP: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3086
GenBank accession number:	BC017021
Gene ID (NCBI):	4223
Full name:	Mesenchyme homeobox 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.