

Labconco® Stackable ScrubberStand™

Maximize your lab space with a Stackable ScrubberStand. By stacking two washers, you'll not only save valuable space, but it will also double your washing capabilities. Experience enhanced productivity without compromising your lab's footprint. The Stackable ScrubberStand is compatible with all Labconco Glassware Washer models.



Features & Benefits

Reduced footprint compared to two undercounter or freestanding units.

Increased washing capacity with two Labconco Glassware Washers at 10.3 ft³.

Strong and stable to support up to 4,000 lbs.


Durable materials such as 1.75" tubular, Glacier White powder-coated steel.

Hardware included to securely mount top washer to stand.


Comparison to Large Capacity Washers

- Increased washing capacity with up to four racks
- Reduced initial starting costs
- Never stop running with redundant washer availability
- Independent application-specific wash programs and cycles

Labconco® Stackable ScrubberStand



Standard ScrubberStand
Catalog Number - 4618000



Vantage ScrubberStand
Catalog Number - 4618001

Glassware Washer Resources



Glassware Washer
Configurator



Glassware Washer
Brochure



Glassware Washer
Information

LABCONCO CORPORATION

8811 Prospect Avenue
Kansas City, MO 64132

(800) 821-5525 | (816) 333-8811

labconco.com



©2024 LABCONCO CORPORATION

Product design subject to
change without notice.

1/11/24