

ELMI V-3 VORTEX



Product Description

VORTEX V-3 is a durable, easily operated, trustworthy addition to every laboratory, and accommodates any test tube up to 50 ml volume.

Features

- 50 to 4,000 RPM
- 4mm rotor vibration amplitude
- Compact design: 5 x 2.4 in.
- 3 Year Warranty
- Amazingly strong
- Hands-free mode
- Simple vortex intensity control
- 2 mixing adapters
- Dome shaped and equipped with anti-vibration system for maximum stability
- Quiet and maintenance-free

Specifications

Rotor Rotation Speed	50 - 4,500
Vortex Rotor Vibration Amplitude	4
Dimensions (D x H)	130 x 58
Power Supply Adapter	12 VDC, 1.5A
Weight	0.7



VRX-02

Mixing Adapter

Total Number of Tubes: 9

Max Tube Capacity (mL): 1.5 - 2ml



elminorthamerica.com

Order: (747) 226-3330

info@elminorthamerica.com

5115 Douglas Fir Rd Suite G | Calabasas, California 91302

© ELMI North America, Inc. August 2016