

# Flexi-Therm Dry Block Heater



A pre-calibrated, accurate sample heater.

## Description

The Flexi-Therm is the next generation of Dry Block Heater.

Used with interchangeable Heated Modules, the Flexi-Therm provides highly accurate sample heating for a range of microplate and tube formats.

Simply load your Module, set a heating temperature and time and wait for the audible alert when your heating program is finished.

## Key Features

- Pre-calibrated Heated Modules
- Accurate sample temperature display
- Compact unit with a small footprint
- Simple user display and touch control interface
- Independent control of sample and lid temperatures to reduce condensation and evaporation
- Memory function for previously used settings

## Accessories

The Flexi-Therm requires at least one Vitl Heated Module which are sold separately.

For the full specifications and range please visit [www.vitlproducts.com](http://www.vitlproducts.com)



## Related Products



**Co-Mix  
V101001**



**Ther-Mix  
V102001**

## Technical Specifications

**Dimensions (mm) W x D x H** 190 x 250 x 80  
(Excluding Heated Module)

**Block heater range** Ambient to 99.9°C  
(Self-limited by some Heated Modules)

**Lid heater range (where fitted)** Ambient to 105°C

**Timer range** 1 minute to 72 hours

**Power Supply** 100 – 240 VAC nominal

## Catalogue number

**V103001**- Flexi-Therm Base Unit

**Vitl**   
LIFE SCIENCE SOLUTIONS  
A member of the ITL group

Please contact [sales@vitlproducts.com](mailto:sales@vitlproducts.com) for more information.

1001\_K\_07\_Iss1

T: +44 1233 638 383 **UK**  
T: +1 804 381 0905 **USA**

**vitlproducts.com**