

IMPROVACUTER® Urinalysis Tubes (MFR470031180)

Features:

1. IMPROVACUTER® Urinalysis Tubes are used for urinary sediment and component analysis.
2. Closed systems that reduce healthcare workers' need to contact with potentially hazardous specimens.
3. The conical Bottom is suitable for routine urine test and microscopic sediment analysis.

IMPROVACUTER® Urinalysis Tubes (MFR471130180)

Features:

1. IMPROVACUTER® Urinalysis Tubes are used for urinary sediment and component analysis.
2. Closed systems that reduce healthcare workers' need to contact with potentially hazardous specimens.
3. The conical Bottom is suitable for routine urine test and microscopic sediment analysis.
4. The preservative tubes contain Sodium Propionate, Chlorhexidine Ethyl p-hydraxybenzoate as additive that provide better sample preservation.

Product Code	Tube Size/Volume	Additive	Color	Closure Type	Label	Material	Pk/Ct Qty
470031180	16x100mm/8ml	/	Yellow	Conventional	Paper	PET (Conical)	50/1000
471130180	16x100mm/8ml	Sodium Propionate, Chlorhexidine Ethyl p-hydraxybenzoate	Red	Conventional	Paper	PET (Conical)	50/1000

Application uses:

Urinalysis Tubes.

