



## Flip-Up Disposable Face Shields DATA SHEET

- 1-in wide plastic headband strap
- Adjustable rear hook & loop fastener fits head circumference 17-in – 26-in
- 14-in W x 8-in H x .020-in thick clear plastic front
- Non-sterile
- Face shield features flip-up and -down with side aluminum rivets



### **PRODUCT:**

Clear Amorphous Polyethylene Terephthalate (APET)

### **COMPOSITION:**

The material is composed of 100% prime grade post-industrial amorphous polyethylene terephthalate.

### **PRODUCT PERFORMANCE STATEMENT:**

Product will perform to Product Specifications provided the following storage and usage conditions are met from the point of receipt.

### **STORAGE CONDITIONS:**

Products are recommended to be kept at temperatures between 40-95°F. Products should not be stored in an area with direct sunlight exposure and should be kept in their original packaging. The storage area should be kept clean and free from contamination or objects that could damage the sheet upon contact. Allow the sheet 24-hours to acclimatize to working conditions before use.

### **SHELF LIFE:**

Six (6) months from receipt of material when stored under the above conditions.

### **USAGE CONDITIONS:**

The recommended usage conditions of thermoformed products produced from sheets 32-140°F.



Properties	Value/Indication	Unit	Tolerance	Test Method
Gauge	10 to 65	mils	+/- 3%	OCTAL Method (QP - 007)
Width	20.0 to 80.0	inches	+/- 0.04	OCTAL Method (QP - 006)
Density (APET)	1.33	g/cm <sup>3</sup>		ASTM D792
Light Transmission	> 86	%		ASTM D1003
Haze	≤ 10.0	%		ASTM D1003
Color L*	> 90	-		OCTAL Method (W04-QC-P03)
Color b*	< 1.0	-		OCTAL Method (W04-QC-P03)
Standard Core ID	6	inches	-	-
Core Alignment	Sheet nom. Width	inches	+ 0/ -0.25	
Roll length (for linear meter orders)	Nominal	Linear meter	Within 10% of nominal length	-
Density Property*	Non-silicone coating w/ anti-fog properties	-		-

**ORDERING INFO:**

Thomas No.	Description	Qty
20A00H653	Face Shield, Flip-Up Disposable, Clear	50/Cs



*Disclaimer: Pricing and offers subject to change without notice.*