


Pluggo™ RHs (Rack Handler-Single) Automated Decapper System for

Cost-effective and Safe Decapping of Vacuum Collection tubes in Hamilton 24 position Racks



- **Bench-top system**
1.5' Long x 3.5' wide x 1.5' Deep
- **No Pneumatics**
- **Connect to any standard outlet,**
Available in 110 VAC 60Hz & in 220VAC 50Hz
- **No more handling of single open tubes**
- **Load a Single Rack straight from the Pluggo to your Instrument**
- **Automatically decaps 24 tubes in 27seconds**



Model 4600 6550

for Hamilton 24 pos. Standard rack (173400) for 16x100mm tubes

Model 4600 6553

for Specialty Hamilton 24 pos. rack (199500) for 16x100mm & 13x100mm tubes & **OPTIONAL 13x75 in LGP Adapters**

Model 4600 6550.75 & 4600 6550.100

available for 13x75 & 13x100 respectively with the use of adapter # 0204449-001 from Hamilton

Protect your lab workers from repetitive motion

PLUGGO IS THE ANSWER!

- **STOP Manually Popping Caps!**
- **Protect Against Repetitive Motion Injuries**
- **Protect Against Aerosols**
- **Improve Laboratory Productivity**