

UBE2R2 Polyclonal Antibody

Catalog number: 14077-1-AP

Size: 37 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

UBE2R2; CDC34B, FLJ20419,

UBC3B, UBE2R2, Ubiquitin

carrier protein R2, Ubiquitin

protein ligase R2



Jurkat cells were subjected to SDS PAGE followed by western blot with 14077-1-AP(UBE2R2 Antibody) at dilution of 1:600

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated UBE2R2 MW:	27 kDa
Observed UBE2R2 MW:	28-30 kDa
Positive WB detected in	Jurkat cells, HeLa cells, mouse spleen tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag5202
GenBank accession number:	BC047584
Gene ID (NCBI):	54926
Full name:	Ubiquitin-conjugating enzyme E2R 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.