

Product Name:	ABT-HC-RFC-15A
Description:	<p>Standard Hydrocarbon Refrigerator &amp; Freezer Combination          ADA Compliant          14 cu. ft. capacity (10 ref./4 freezer)          Two exterior solid swing doors, right hinged, reversible          Mechanical thermostat          Refrigerator temperature range: 0°C to 10°C          Freezer chamber: -15°C to -10°C with two auto defrost cycles          warming to 0°C per day          Two adjustable glass shelves (refrigerator) and two-gallon door bins,          plus one half gallon door bin          1 full width glass shelf (freezer) and 2 fixed door bins          Magnetic door gasket for positive seal          Integrated door handles          Leveling legs          Interior LED refrigerator light          Two probe access ports (3/8"); one, refrigerator/one, freezer          White textured cabinet and door          High density urethane foam cabinet and door insulation          HFC-free refrigerant (R600a), foam insulation and packaging          Two year parts and labor warranty, plus an additional three year          compressor parts warranty          Exterior dimensions: 27 3/8"W x 29"D x 60 3/8"H          Auto defrost          115V, 60 Hz, 1.1 Amps          Shipping weight: 135 lbs.          DISCLAIMER: Auto defrost freezers incorporate an electric heater          installed on the evaporator that assist in the removal of the frost and          ice. Although the heater helps speed up the frost and ice removal, it          can raise the temperatures within the freezer.</p>