

Product Name:	AJ-IV CS
Description:	<p>ABS Auto Jr. with Fiberboard Package System Total Vial Capacity (2 mL vials): 6000 Total Box Capacity: 60 Rack Capacity: 6 Boxes per rack: 10 Vial per box (2 mL): 81/100 LN2 Capacity (Liters): 165 Static Evaporation Rate (Liter/Day): 0.84 Static Holding Time (days): 194 Power requirement: 24 Vac Neck Diameter (mm): 215.9 Exterior Height (mm): 990.6 Exterior Diameter (mm): 683.3 iSentry Controller</p> <p>Controller features:</p> <ul style="list-style-type: none"> ○ Temperature display (Celsius) ○ Level accuracy with 4 direct read sensors ○ Large format LED display ○ Alarm ring back ○ Redundant level sensor system ○ Remote monitoring system available ○ Liquid level indicator ○ Alarm package (audible, visual, and remote) including level (high and low), high temperature, and sensor failure <p>UL-Listed transformer Stainless Steel armor-plated transfer hose State of the art hanging rack design Cast aluminum roller base Lid with locking provision & protective boot Fiberglass restricted opening neck design One-year warranty on all parts and labor Three-year vacuum warranty</p> <p>Package System Includes:</p> <ul style="list-style-type: none"> - (60) 81 or 100 Cell Fiberboard Boxes & Dividers - (6) 10-2 Racks - (1) Pair CryoGloves - (1) LN₂ Measuring Ruler