

## OSBPL3 Polyclonal Antibody

Catalog number: 12417-1-AP


Size: 30 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

OSBPL3; OSBPL3, ORP3, ORP  
3, OSBP3, KIAA0704



Y79 cells were subjected to SDS  
PAGE followed by western blot with  
12417-1-AP(OSBPL3 antibody) at  
dilution of 1:400

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Calculated OSBPL3 MW:	887aa,101 kDa
Observed OSBPL3 MW:	101kd
Positive WB detected in	Y79 cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3044
GenBank accession number:	BC017731
Gene ID (NCBI):	26031
Full name:	Oxysterol binding protein-like 3

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.