

# Equipment & Instrumentation

Fall / Winter 2023

## Source Guide



**Thomas**  
Scientific

## Featuring Products & Solutions

[ThomasSci.com](http://ThomasSci.com)

# Welcome!

Thomas Scientific is a distributor of laboratory supplies, equipment and chemicals. The products in this catalog represent only a small offering from our extensive catalog and website. We're here to make your purchasing process as simple as possible. For additional models, sizes, packaging, and accessories not shown in our catalogs or if you just need help selecting a product, contact us at 833.544.SHIP (7447) or visit us on the web at ThomasSci.com.



Thomas Scientific accepts VISA, Master Card, American Express and the IMPAC card for government buyers. POs accepted for registered accounts.

## Our Mission

To help science create a better world by providing quality products with personal, timely and efficient service.

Delivering Quality Scientific Supplies  
Worldwide Since 1900

Special offers and pricing are valid through December 31, 2023, unless otherwise noted.

**NOTE:** Offers in this catalog are void where prohibited by law or company policy. Shipping charges are not included. Prices and offers are subject to change without notice and are not to be combined with other offers/discounts or contract pricing. Valid in USA only. Dealer/reseller sales not applicable. For web orders, enter the promo code, if applicable, at time of checkout. For phone orders, please supply the promo code to the customer service representative when ordering. Discount pricing shown must be specified when placing order. Returned products will be subject to billing for Free Items at list price. Free items in 'Buy/Get 1 Free' offers will be the item of lowest value unless otherwise noted. Giveaways and free products may be shipped separately. Allow 4-6 weeks for delivery.

To view more of Thomas Scientific's latest publications and promotions, visit [ThomasSci.com/deals](https://www.thomassci.com/deals) Resources Tab

or Click Here ►►



## Brands Featured In This Issue

**Benchmark**  **ECOLAB** 

 **BRANDTECH** 

  
Horizon Scientific, Inc.



**labForce**   
powered by Thomas Scientific

**LBI** Laboratory Builders, Inc.

  
**METTLER TOLEDO**

**SARTORIUS**

**SCHULER SCIENTIFIC**

**thermo**  
scientific



**Thomas Scientific**  **PODCAST**

Listen While You Work, by checking out Thomas Scientific's latest Podcast.

Our podcasts deliver high quality, valuable content to help streamline your applications and are available on all of your favorite platforms.


   

Visit our homepage for our latest podcast or visit:  
[ThomasSci.com/podcast](https://www.thomassci.com/podcast)





The **RC 6** chemistry-HYBRID Vacuum Pump combines rotary vane and diaphragm vacuum pump technology to support applications in the fine vacuum range with dramatically reduced service demands. In the **RC 6**, a chemistry-design diaphragm pump is used to back the rotary vane pump.



[Click Here](#) ►



The Labconco Logic+ Class II A2 Biosafety Cabinet excels in user comfort and functionality. Industry leading advances in ergonomic design have shaped the Inclination™ Technology Suite that includes a 10° angled sash, patented curved air grille with built in armrest, fully accessible controls, and eye-level mounted color display with MyLogic™ operating system.



◀ [Click Here](#)

## Go Beyond Routine Weighing with XPR and XSR Excellence Balances

Samples are the core of every analysis. Achieving reliable results depends on the accuracy of the initial sample weighing, and METTLER TOLEDO XPR and XSR analytical and precision balances are ready to meet the challenge of handling rare, expensive or demanding samples in complex applications.

Our balance portfolio is the answer to your everyday weighing challenges and enables you to achieve uncompromised weighing excellence. Learn more today by visiting [www.thomassci.com/mettlertoledo](http://www.thomassci.com/mettlertoledo).



[Click Here](#) ►





## SmartReader™ UV-Vis Microplate Reader



The Accuris SmartReader™ UV-Vis (model MR9610) is a welcome addition to any lab that is routinely measuring concentration or absorbance in 96 and/or 384 well plates. Through absorbance readings from both the ultraviolet and visible spectrum range, it can determine the concentration of samples, conduct

positive or negative detection of antibodies and antigens, and perform nucleic acid/protein quantification in samples. The instrument's 10-inch touch screen, intuitive software and graphical interface makes it easy to use as a standalone instrument (a separate computer and additional software is not required).

Visit [ThomasSci.com](http://ThomasSci.com) Keyword **22A00A218** or [Click Here](#) ▶▶



## Automated and Adaptable

Discover Ecolab's range of bio-decontamination systems using Bioquell Hydrogen Peroxide Sterilant\*



Visit [ThomasSci.com](http://ThomasSci.com) and search for Bioquell



GLOBAL EXPERTS IN CONTAMINATION CONTROL

© 2023 Ecolab USA Inc. All rights reserved. 28APR23/NA



\*Bioquell Hydrogen Peroxide Sterilant (EPA Registration Number: 72372-1-86703) Use Bioquell products safely. Always read the label and product information before use.



# Perform convenient field separations!

Compact, lightweight, and designed for easy transport.

The Globe Scientific GCC-P **Portable Clinical Centrifuge** utilizes an included 12VDC car power adapter and has all the features required to perform routine clinical applications in the field. A weighted base ensures the unit remains stable in field environments.

- Suitable for numerous clinical and research applications
- Operates using an included 12VDC car adaptor
- Fast acceleration and two deceleration speeds
- Speed range of 300-5000rpm
- Independent control of speed or RCF & time
- Audible alert function
- LCD display and automated lid lock



Accommodates 1.5mL to 5mL, 7mL, 10mL and 15mL tubes.



**globe**<sup>™</sup>  
SCIENTIFIC



Click Here



21A00D391

**COREPOINT**  
SCIENTIFIC

Lab Freezers that  
**Save Energy**

Corepoint<sup>®</sup> Scientific laboratory freezers use hydrocarbon refrigerants for lower energy consumption and a cleaner planet.



Visit [ThomasSci.com](http://ThomasSci.com)  
Search 21A00D391 or Click Here




**Thomas Quick Ship Program:**

**The Complete Source For**

**Lab Furniture & Accessories You Need Now.**

***Click [HERE](#) to view our  
online Quick Ship Catalog***

**Thomas  
Scientific**



Quick Delivery Lab Furnishings Catalog

**INSTOCK**  
Quick Delivery Program

**Thomas Scientific's laboratory furniture team offers:**

- A quick delivery furniture program with complete labs shipping within 14 business days\*
- An extensive inventory of high quality furnishings & components
- Clear and easy project communication

**Services include:**

- Design, project management and factory trained installation

**Thomas  
Scientific**

# SARTORIUS Sartolab® Vacuum Filtration Units



## Special Offer!

Save Up to 35% on Sartolab® Vacuum Filtration Units

- Highest flow rates with a large filtration area
- No loss of protein with a low protein affinity membrane
- Low hold-up volume due to an optimized membrane support
- Can be used as stand-alone system or with the Sartolab® Multistation, which allows parallel filtration of up to 6 samples

Special offer valid only for the US until December 31st, 2023 and cannot be combined with other offers.

[Click Here](#) ▶▶



## labForce

powered by Thomas Scientific

### TS-T Series Laboratory Balances

- Readability options ranging from 0.01 mg to 0.01 g
- Internal Calibration
- Option to purchase with an ISO 17025:2017 Calibration Certificate included
- 5" color touchscreen
- Communication interfaces including 2xRS232, USB-A, USB-B, Ethernet, Wi-Fi®



◀◀ [Click Here](#)



**thermo**  
scientific

# Sorvall™ BIOS Centrifuge Series

Centrifuges designed for reliable,  
high-throughput bioprocessing

- Outstanding capacity to maximize productivity
- High performance for application flexibility
- Enhanced ergonomics with push-button convenience
- Simplified operation and space-saving design



Visit [ThomasSci.com](http://ThomasSci.com)  
Search "Sorvall" or [Click Here](#) ▶▶



## thomas experts series

delivers the experts to you

**Ask your Thomas Rep about the Thomas Experts Series!**

Every month the **Thomas Experts Series** will bring emerging trends, new applications, and workflow solutions directly to you, wherever you are. Our Thomas Talk Experts will discuss a variety of timely topics, ranging from Chemical & Laboratory Safety to Automation and Environmental Science Technologies.

Our rich virtual environment is filled with PDFs, web links, videos, incentives, surveys, and of course our Thomas Expert Talks! Clicking, downloading, filling your virtual briefcase, registering, and attending a Thomas Talk enters you into multiple drawings, including the Thomas Experts Series monthly Beats Wireless Noise Cancelling Bluetooth Earbuds.



### Expert Items:

- Expert welcome greeting
- Zoom, drawings, survey links
- Literature pdfs, web links, videos a resource rich information kiosk
- Monthly and up coming topics and Thomas Experts Series Schedule
- Thomas Specials, Services, Requests

**Scan to register for 1  
or any number of our  
24 Thomas Expert  
Talks or [Click Here!](#)**

Connect With Us:



**Thomas**  
Scientific