

PCTAIRE1 Polyclonal Antibody

Catalog number: 10102-1-AP

Size: 35 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:


PCTK1; CDK16, FLJ16665,

PCTAIRE, PCTAIRE motif

protein kinase 1, PCTAIRE

protein kinase 1, PCTAIRE1,

PCTGAIRE, PCTK1



mouse brain tissue were subjected to SDS PAGE followed by western blot with 10102-1-AP(PCTK1 antibody) at dilution of 1:600

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated PCTAIRE1 MW:	56 kDa
Observed PCTAIRE1 MW:	56-63kd
Positive WB detected in	Mouse brain tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag0151
GenBank accession number:	BC001048
Gene ID (NCBI):	5127
Full name:	PCTAIRE protein kinase 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.