

Product Datasheet

Product code: **8WF.SA.1**

Product type: Wafer Tweezers

Made in: Switzerland

Description: Type 8WF wafer tweezers are suitable to delicate applications which require holding a 6-inch flat component firmly. The line serration provides the handles with extra grip. The bent paddle helps to hold the component still during positioning.

Tip description: handling 6", top fingers, stepped bottom paddle


Body material: Anti-Magnetic Anti-Acid Stainless Steel

Picture:



Tweezer dimensions:	OAL	135.0 mm – 5.315 inch
	Width	28.5 mm – 1.122 inch
	Height	12.0 mm – 0.472 inch

Tip details:



Tip dimensions (Wafer tweezers):	A (tip width)	28.5 mm – 1.122 inch
	B	12.0 mm – 0.470 inch
	C	2.2 mm – 0.080 inch
	D	3.5 mm – 0.140 inch
	E	2.5 mm – 0.100 inch