



FEEL THE SW DIFFERENCE[®]

DESIGNED TO ENHANCE WORKER PERFORMANCE.

PowerForm[®]

PF-95BL

POWDER-FREE ECOTEK SUSTAINABLE
NITRILE EXAMINATION GLOVES WITH DRY CORE*

KEEP HANDS COOLER AND MORE COMFORTABLE WITH PREMIUM EXAM QUALITY NITRILE GLOVES, FEATURING UNIQUE DRY CORE MOISTURE MANAGEMENT TECHNOLOGY.

- Dry Core[®] Moisture Management Technology pulls sweat away from the skin
- EcoTek[®] Sustainable Technology*
- Silicone-free for critical environments

ANSI Level 3
Abrasion
Protection

Fully textured
surface for wet
and dry grip

Tear-resistant
beaded cuff for
easy donning



PF **PF-95BL**

WEIGHT	LENGTH	THICKNESS	FEATURES	COLOR
 8.5g (large size)	 9.5"	 6.1 mil	Dry Core [®] EcoTek [®]	 BLUE

Ecotek



SW[®]
33278 Central Avenue, Unit 102, Union City, CA 94587, USA
Tel: +1.510.429.8692 | Fax: +1.510.487.5347

Trademarks and registered trademarks are the property of SW and its affiliates.
©2023 SW. All rights reserved.



TECHNICAL DATA



PF-95BL POWDER-FREE ECOTEK SUSTAINABLE NITRILE EXAMINATION GLOVES WITH DRY CORE*



PRODUCT DESCRIPTION

SPECIFICATIONS

PRODUCT	PowerForm®		
GRADE	Powder-Free Exam		
MATERIAL	Nitrile		
THICKNESS	FINGER	6.9mil	0.18mm
	PALM	6.1mil	0.16mm
	CUFF	5.1mil	0.13mm
COLOR	OUTER	Blue	
	INNER	White	
SURFACE	OUTER	Fully Textured	
	INNER	Dry Core®	
LENGTH	9.5"	242mm	
WEIGHT	8.5g		
FIT TYPE	Ambidextrous		
CUFF TYPE	Beaded		
WATERTIGHT AQL	1.5		
ORIGIN COUNTRY	China		

PHYSICAL PROPERTIES

	PRE-AGING	POST-AGING
TENSILE STRENGTH:	20 MPa	19 MPa
ELONGATION:	560%	550%

INTEGRATED TECHNOLOGIES

DRY CORE Moisture Management Technology: Absorbent lining wicks away moisture for added comfort in extended use.

ECOTEK Sustainable Technology: Product has been tested per ASTM D5526*



* Specifications based on size Large glove

STANDARDS & CERTIFICATIONS

EN ISO 21420:2020
EN 455-2000-2015
ASTM D6319-19
ASTM F1671-13
ASTM F739-12
Reach EU 1907/2006
EU 2016/425

ANSI/ISEA 105-2016
CAT III Complex Design
FDA 21 CFR 177.2600



Please visit swwglobal.com/symbols for information on standards used on this document.

ORDERING

ITEM NO.	SIZE	PALM WIDTH	PACKAGING
N106551	XS 5½ – 6	3.0" 76mm	50 gloves per box 10 boxes per case
N106552	S 6½ – 7	3.4" 87mm	50 gloves per box 10 boxes per case
N106553	M 7½ – 8	3.8" 97mm	50 gloves per box 10 boxes per case
N106554	L 8½ – 9	4.3" 109mm	50 gloves per box 10 boxes per case
N106555	XL 9½ – 10	4.5" 114mm	50 gloves per box 10 boxes per case
N106556	2XL 10½ – 11	4.9" 125mm	50 gloves per box 10 boxes per case

APPLICATIONS

- For Medical Use
- Chemical Handling
- Clinical
- Laboratory
- Liquid Handling
- Pharma Production

Keep out of sunlight. Store in a cool, dry place. Keep away from sources of ozone or ignition.



R_N10655X0323



SW®
33278 Central Avenue, Unit 102, Union City, CA 94587, USA
Tel: +1.510.429.8692 | Fax: +1.510.487.5347

Trademarks and registered trademarks are the property of SW and its affiliates. ©2023 SW. All rights reserved.



*EcoTek gloves tested to ASTM D5526-12 Standard Test Method for Determining Anaerobic Biodegradation of Plastic Materials Under Accelerated Landfill Conditions. Results are 92.6% biodegradation in 60% solids landfill at 945 days. Future results cannot be predicted/extrapolated. If you need more information on the test reports, send email to cs@swwglobal.com.