

RPA3 Polyclonal Antibody

Catalog number: 10692-1-AP

Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

RPA3; REPA3, Replication

factor A protein 3, replication

protein A3, 14kDa, RF A

protein 3, RP A p14, RPA14,

RPA3

Applications

Tested applications:	ELISA
Species specificity:	Human; other species not tested.
Calculated RPA3 MW:	13 kDa
Observed RPA3 MW:	

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag1065
GenBank accession number:	BC005264
Gene ID (NCBI):	6119
Full name:	Replication protein A3, 14kDa

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.