

# Lab Specimen Transport Bags

## Lab Guard® Lab Specimen Transport Bags — Zip Closure, Utility (No Biohazard Print)

Lab Guard® UTILITY Specimen Bags are not printed with a biohazard symbol for hospitals and labs with in-facility processing and no transportation outside of the facility.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging	Qty/cs	
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
SBL2X69	49-90	6x6	Single Zip	Clear	NONE	NONE	NONE		w/flap	Zip Bag	10/100	1000
		6x9	Double Zip	Clear/Blk & Org	NONE	NONE	NONE	TearZone®	no flap	DispBag	20/50	1000
	49-96	6x9	Single Zip	Clear	NONE	NONE	NONE		w/flap	Zip Bag	10/100	1000
	49-93	8x8	Single Zip	Clear	NONE	NONE	NONE		w/flap	Zip Bag	10/100	1000



49-90



SBL2X69

## Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure

Speci-Gard® Specimen Bags feature the innovative, easy-to-use, liquid-tight PRESS AND CLOSE® permanent adhesive closure, TearZone® EZ open system, DispenserBag® convenience packaging, and in some cases, a cost-saving destroyable biohazard symbol.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging	Qty/cs	
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
95-CMS59		5x9	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
UF95-66SPBIO		6x6	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	no pouch	DispBag	100	2000
UF95-600BIO		6x10	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
UF95-600BPM		6x10	Adhesive	Clear/Blk & Red		X	Standard	TearZone®	w/flap	DispBag	100	1000
95-CMS810		8x10	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
UF95-610BIO		10x10	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
UF95-1318BIOA		13x18	Adhesive	Clear/Blk & Red		X	Standard	TearZone®	w/flap	Zip Bag	50	500
UF95-1517BIO		15x17	Adhesive	Clear/Blk & Red		X	Standard	TearZone®	w/flap	Zip Bag	50	500



UF95-600BIO



UF95-600BPM



SPEC-PNYRT



UF95-600STB



UF95-600STBR

# Lab Specimen Transport Bags

## Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure - STAT

Speci-Gard® STAT Specimen Bags are available in red and yellow film colors with all featuring the easy-to-use, liquid-tight, PRESS AND CLOSE® permanent adhesive closure, TearZone® EZ open system, cost saving destroyable biohazard symbol and DispenserBag® convenience packaging.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging	Qty/cs	
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
SPEC-PNYRT		6x10	Adhesive	Red Tint/Blk		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
UF95-600STBR		6x10	Adhesive	Red Opaque/Blk		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
UF95-600STB		6x10	Adhesive	Yel Opaque/Blk		X	Destroyable	TearZone®	w/flap	DispBag	100	1000

## Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure & Absorbent Pad Insert (3x5 Pad)

Speci-Gard® Specimen Bags with absorbent pads are available with either single adhesive pouch (UF95-BDES) or DUAL-adhesive pouch (SPECI-CC). Item SPECI-CC secures both the specimen and the paperwork in separate pouches so that neither can be tampered with, without destroying the bag (chain-of-custody).

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging	Qty/cs	
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
UF95-BDES		6x10	Adhesive	Clear/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
SPECI-CC		6x10	Adhesive	Clear/Blk & Red		X	Standard	TearZone®	Adhesive	DispBag	100	1000



UF95-BDES



SPEC-PNYB



SPEC-PNYG



SPEC-PNYP

## Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure - Color Codes

Speci-Gard® Color-Code PRESS AND CLOSE® Adhesive Closure Specimen Bags feature a range of film colors to help hospitals and labs segregate specimens with specific processing needs.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging	Qty/cs	
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, Rm Temp						
SPEC-PNYB		6x10	Adhesive	Blu Tint/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
UF95-600BBL		6x10	Adhesive	Blu Opaq/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
SPEC-PNYG		6x10	Adhesive	Grn Tint/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
SPEC-PNYP		6x10	Adhesive	Pur Tint/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000
SPEC-PNYY		6x10	Adhesive	Yel Tint/Blk & Red		X	Destroyable	TearZone®	w/flap	DispBag	100	1000

SPEC-PNYY



UF95-600BBL



UF95-600



# Lab Specimen Transport Bags

## Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure, Utility (No Biohazard Print)

Speci-Gard® UTILITY Specimen Bags are not printed with a biohazard symbol for hospitals and labs with in-facility processing and no transportation outside of the facility. All bags offer TearZone® EZ-Open Feature and DispenserBag® Convenience Packaging.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging	Qty/cs	
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, RmTemp						
UF95-66SP		6x6	Adhesive	Clear/Blk & Red	NONE	NONE	NONE	TearZone®	w/flap	DispBag	100	1000
UF95-600		6x10	Adhesive	Clear/Blk & Red	NONE	NONE	NONE	TearZone®	w/flap	DispBag	100	1000
UF95-610		10x10	Adhesive	Clear/Blk & Red	NONE	NONE	NONE	TearZone®	w/flap	DispBag	100	1000

## Speci-Gard® Lab Specimen Transport Bags: Permanent Adhesive Closure, UV Protection

Speci-Gard® style PRESS AND CLOSE® Adhesive Closure Bag features opaque gray white film that prevents light from entering the bag, protecting light sensitive contents. Access to contents is available with EZ-Tear Perforation.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging	Qty/cs	
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, RmTemp						
SPECI-UV		6x10	Adhesive	Wht Opq/Blk & Red	NONE	NONE	Standard	Perforation	w/flap	DispBag	100	1000

## Lab Specimen Transport Bags: Permanent Adhesive Closure, Foldover Closure

Lab specimen bags feature a traditional foldover permanent adhesive closure and score line for opening during processing.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging	Qty/cs	
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, RmTemp						
	59-89	6x9	Adhesive	Clear/Blk & Red	NONE	NONE	Standard	Score	no flap	Zip Bag	100	1000
	4060	6x10	Adhesive	Clear/Blk & Red	NONE	NONE	Standard	Score	no flap	Zip Bag	100	1000
	59-93	10x10	Adhesive	Clear/Blk & Red	NONE	NONE	Standard	Score	no flap	Zip Bag	100	1000

## Lab Specimen Transport Bags: Adhesive Closure, Specialty Applications

Item 1130 features rugged welded side seams and permanent adhesive closure for pneumatic tube systems.

Item 39-97 features a reclosable adhesive closure for worry-free post-closure access.

Minigrip Code#	Medegen Code#	Size	Closure	Bag/Print Color	Selector Box Print		Biohazard Symbol	Opening	Doc Pouch	Packaging	Qty/cs	
					Stat, Frozen, Refrig, RmTemp	Frozen, Refrig, RmTemp						
	1130	6x9	Adhesive	Clear/Blk & Red	NONE	NONE	Standard	Tear Notch	no flap	Bag	100	1000
	39-97	6x10	Reclosable Adh	Clear/Blk & Red	NONE	NONE	Standard	Pull Open	no flap	Zip Bag	100	1000
	39-99	8.5x11	Reclosable Adh	Clear/Blk & Red	NONE	NONE	Standard	Pull Open	no flap	Zip Bag	100	1000



SPECI-UV



4060



1130



39-97