


UltraDose® Germicidal Cleaner

Features:

- E.P.A. registered Product
- For use in an ultrasonic cleaner
- Disinfects instruments and the ultrasonic bath
- Cleans non-porous instrument and appliance surfaces, apparatus and equipment
- Performs as an ultrasonic cleaning solution, disinfectant, deodorant, germicidal, fungicidal, virucidal and tuberculocidal all-in-one
- Concentrated; one ounce makes one gallon.
- Non-irritating and fragrance free



Thomas No.	Description	Quantity
20A00H588	UltraDose® Germicidal cleaner - 16oz bottles	6/Case

Disclaimer: Pricing and offers subject to change without notice.