

## MYO7B Polyclonal Antibody

Catalog number: 14467-1-AP

Size: 41 µg/150 µl


Source: Rabbit

Isotype: IgG


Synonyms:

MYO7B; DKFZp686A08248,

MYO7B, myosin VII B



Immunohistochemical of paraffin-embedded human colon cancer using 14467-1-AP (MYO7B antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemical of paraffin-embedded human colon cancer using 14467-1-AP (MYO7B antibody) at dilution of 1:50 (under 10x lens)

### Applications

Tested applications:	ELISA, IHC
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated MYO7B MW:	242 kDa
Observed MYO7B MW:	
Positive IHC detected in:	Human colon cancer tissue
Recommended dilution:	IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag5829
GenBank accession number:	BC035615
Gene ID (NCBI):	4648
Full name:	Myosin VII B

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.